plusmed

health management

Model: **pM-BS 01**

ELECTRONIC BABY SCALE
ELEKTRONİK BEBEK BASKÜLÜ

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EN Instruction Manual

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Batteries – 4 x AA 1.5V alkaline (Included)

Accuracy: ± (3%+5g)

Working temperature: 10-40 °C/ Working temperature: 0-104°F

Storage temperature: -10-50 °C / Storage temperature: 41-122°F

Humidity: <=90% relative humidity

Advice

- If the scale is not used frequently, remove the batteries.
- Place the scale on a flat, even surface.
- Do not overload the scale.
- Avoid knocking or dropping the scale.
- When cleaning the scale do not use harsh chemical or rough cleaning liquids.
- Do not immerse in water or any other liquid.
- If the scale functions irregularly, remove the batteries, wait for a few minutes and replace the batteries.

Weighing

- Insert 4 new batteries in the battery compartment under the scale. Replace the battery cover.
- Put the scales on a flat, even surface. A small blanket can be laid on the weighing tray.
- Press the **ON/OFF** key to start.
- The display shows **88888**.
- Wait for the display to show **0.000** kg
- Place the baby at the centre of the tray.
- The weight is shown for 60 seconds.
- Press the **ON/OFF** button to turn off the scale manually. The scale turns off automatically after 60 seconds.

Tare – Add-and-weigh

While the scale is turned on and at **0.000** kg, if you place an object on the scale, e.g. blanket. The blanket weight is shown. Press the **TARE** key to return to **0.000** kg. You can now place the baby on the scale and weigh.

Minus sign - While the scale is turned on and at **0.000** kg, if you place an object on the scale, e.g. blanket. If you immediately remove the blanket, a minus sign will show in the display. Press the **TARE** key to come back to **0.000**kg.

Recall previous measurements

While the scale is turned on, press the **SET** key. The previous measured weight will be displayed. Continue to

press the **SET** key and up to 10 previous readings will be displayed. The weight shows for 3 seconds and then returns to weighing mode.

To change the unit of measurement for weight:

While the scale is turned on, press the **UNIT** key to change the unit of measurement. The factory setting is kg and °C.

To change the unit of measurement for temperature:

Press the **ON/OFF** key to start.

While the display shows 88888 press the unit key. This changes the temperature only.

LCD displays

Lo - If "Lo" shows in the display, the batteries are low. Replace with new batteries.

0-Ld - The scale is overloaded. The weight is more than 20.10kg. Remove the weight immediately.

UnSt - If you try to return the scale to **0.000**kg and the scale is unstable, **UnSt** will show in the display. Make sure the scale is on a flat, even surface and start again.

Length measurement

There is a flexible tape measure built into the right side of the scale. Pull this out to measure the length of the baby. After using, press the plastic button to stow away the measuring tape.

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