



## pM-T01K | Infrared Ear Thermometer

This ear thermometer has been clinically tested in several big teaching hospitals based on the protocol of EN 12470-5:2003 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading. The continuous ear scanning mode views the entire ear canal and looks for the highest temperature that can be detected. Simply press and hold the start button for 3 seconds, gently rotate the probe inside the ear canal in a circular manner and the most accurate reading will appear.

- 1 1 second measurement or continuous scanning mode
- 2 Probe Cover free
- 3 Fever alarm
- 4 9 memories
- 5 °C / °F switchable
- 6 Backlight
- 7 Auto shut-off
- 8 Low battery indicator

